

Eisner™

A division of Multigate
Medical Products Pty Ltd



Scissor Inserts

Eisner™ laparoscopic
instruments deliver
unmatched durability,
comfort and efficiency
to support a seamless
performance every time.

1993

POWERING
LAPAROSCOPIC
SURGERY
SINCE 1993



DESIGN

Our instruments are designed for optimal balance, feel, comfort and responsiveness to deliver a seamless connection between the hand and the surgical site.

SAFETY & RELIABILITY

The craftsmanship and the quality of materials we select combine to deliver high performance, electro-safe instruments that deliver unrivalled durability, longevity and reliability.



• **Easy to re-process** Our instruments are designed to simplify re-processing.

• **Clinical performance** Our hand-held instruments are modular allowing you to build cost effective instrument trays whilst maximising patient safety and clinical performance.



INSERTS - SCISSORS

We offer a range of straight and curved blade scissor inserts.



Inserts - Scissors (5mm)

Scissors with Insulated Tube	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Metz Curved Scissor (Tri-Cut) with Insulated Tube & L'nut	300-600-177	300-600-177.45	15mm Blade	Double Action Medium Curve Blade Sheath and Locking Nut Included	1

Curved Blade	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Curved Multicut Scissors	MC2	—	11mm Blade	Double Action Mild Curve Maximum of 10 Procedure Use Recommended	1
	5mm Curved Multicut Scissors (Pack of 5)	MC2000	—	11mm Blade	Double Action Mild Curve Blade	1 Pack

	5mm Metz Strong Curve Baby Scissors Insert	552030	—	14mm Blade	Double Action Strong Curve Metzenbaum Blade	1
	5mm Metz Strong Curve Scissors Insert	552112	—	16mm Blade	Double Action Strong Curve Metzenbaum Blade	1
	5mm Metz Medium Curve Scissor Insert	400067	401067	12mm Blade	Double Action Medium Curve Metzenbaum Blade	1
	5mm Metz Strong Curve NPHM Scissor Insert	400072	401072	15mm Blade	Double Action Non-protruding Hinge Mechanism	1
	5mm Micro Curved Scissors Insert	400572	401572	7mm Blade	Double Action Short Petite Curved Blade	1

INSERTS - SCISSORS (5MM)

Straight Blade	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Metz Straight Scissor Insert	400069	401069	12mm Blade	Double Action Straight Metzenbaum Blade	1
	5mm Micro Straight Scissors Insert	400570	401570	7mm Blade	Single Action Short Blade Wide Opening Feature	1
	5mm Metz Straight Scissors Insert	400068	401068	15mm Blade	Double Action Straight Longer Metzenbaum Blade	1
	5mm Straight Multicut Scissors	MC1	—	11mm Blade	Double Action Straight Blade	1
	5mm Straight Multicut Scissors (Pack of 5)	MC1000	—	11mm Blade	Double Action Straight Blade	1 Pack

Other Scissors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Hook Scissor Insert	400060	401060	8mm Blade	Single Action Hook Style Blade	1
	5mm Knot Pusher Scissor Insert	400331	401331	16mm Blade	Double Action Knot Pusher Tip & Scissor in Unique Molding	1

SUPPORT

With Eisner™ you don't get a call centre, just a dedicated team of specialists ready to tailor your Eisner™ experience.

Experience the difference

Book an obligation free trial

Our team can organise an obligation free, trial assessment of your instruments so you can see and feel the difference before you buy.

Trade-in / replacement program

Manage your tray upgrades

If you are considering decommissioning an instrument our trade-in/replacement service will assist you to address one off replacements or the gradual upgrade of instrument sets. We accept Eisner™ and non-Eisner™ instruments.

Care

Optimise Your Tray Performance

When you choose Eisner™ products our team are available to provide free comprehensive in-service training, support and tray assessments to guide servicing and repair.





Enhance your every move

Warranty | Eisner™ offers a 12 month warranty against manufacturing defect on all metal components of the equipment. Damage and / or breakages caused by improper use or abuse of the equipment are not covered by warranty. Normal wear and tear is not covered. Warranty is null and void where components are mixed with non-Eisner™ instruments.

Service & Repairs | Instruments returned to us for repair/replacement will be assessed and a written quotation provided for you to make a decision on further action. In most cases for instruments beyond repair we can offer a trade-in replacement. Eisner™ keeps factory original spare parts in stock. All repairs should be returned to Eisner™ for repair in order to preserve the integrity of your instruments.

Recommended Lifespan of Instruments | Whilst Eisner™ laparoscopic instruments will function for many years, we recommend replacement prior to any obvious signs of wear and stress which could lead to an issue during surgery. We recommend a lifespan of 3 years on our handles and inserts and suggest that after that time, they are traded in on a new component. Insulated Tubes: By virtue of the fact that insulated tubes will be replaced once the insulation is damaged or compromised in any way, it is not envisioned that the lifespan of a sheath will be beyond 12 months. Instruments should be taken out of circulation immediately if it is found that identifying labelling is no longer visible.

EisnerTM

A division of Multigate
Medical Products Pty Ltd

Pacifica Business Park
Unit 18, 97 Canterbury Road
Kilsyth VIC 3137
Australia
T +61 (0)3 9728 4770
F +61 (0)3 9728 6770
eisner@multigate.com.au
www.eisner.com.au

For an obligation
free trial or for more
information call
us today on:
(03) 9728 4770
eisner.com.au