


Klein® Bulldog Clamps and Applier

Bulldog clamps are used to securely grip and hold structures, tissue or vessels, stop blood flow to an organ or prevent excessive bleeding in an area. The Eisner range contains 10 different styles of bowel, vascular and trocar clamps so you can confidently select an ideal clamp for the task at hand.


CLAMPS KEY FEATURES AND BENEFITS

- ✓ Easy application
- ✓ Rapid clamp removal
- ✓ Debakey atraumatic serrations
- ✓ Gold dipped to improve intra-operative visibility
- ✓ Easy to clean. No spring or intrinsic working mechanism



APPLIER KEY FEATURES AND BENEFITS

- ✓ 3 Piece dismantlable, rotating, Ø10mm, 33cm WL stainless steel applier
- ✓ Quick and easy assembly/disassembly for cleaning and sterilisation
- ✓ Spring action handle provides automatic clamp closure and applicator disconnection
- ✓ Holding bar holds the clamp on the applier while positioning
- ✓ AS/NZS 4187 compliant



Bulldog Clamp - Klein®	DESCRIPTION	RE-ORDER 33cm	SIZE	FEATURES
	Endoscopic Klein® Bulldog Clamp Applier	300-410-775	Ø10mm	Rotating, 33cm length Applier for Detachable Clamps
	Endoscopic Klein® Bulldog Bowel Straight Clamp	300-100-771	45mm Jaw	Straight Mid-length Clamp
	Endoscopic Klein® Bulldog Bowel Straight Clamp	300-100-772	60mm Jaw	Straight Long-length Clamp
	Endoscopic Klein® Bulldog Vascular Straight Clamp	300-100-773	17mm Jaw	Straight Short-length Clamp
	Endoscopic Klein® Bulldog Bowel Straight Clamp	300-100-774	25mm Jaw	Straight Mid-length Clamp
	Endoscopic Klein® Bulldog Vascular Curved Clamp	300-100-776	17mm Jaw	Mild Curve Short-length Clamp
	Endoscopic Klein® Bulldog Vascular Curved Clamp	300-100-777	25mm Jaw	Mild Curve Mid-length Clamp
	Endoscopic Klein® Bulldog Bowel Curved Clamp	300-100-778	45mm Jaw	Mild Curve Long-length Clamp
	Endoscopic Klein® Bulldog Bowel Curved Clamp	300-100-779	60mm Jaw	Mild Curve Extra Long-length Clamp
	Endoscopic Klein® Bulldog Satinsky Back Curved Trocar Clamp	300-100-786	60mm Jaw	Long Back Curved Clamp To be use with a Ø12.5mm Trocar System
	Endoscopic Klein® Bulldog Satinsky Forward Curved Trocar Clamp	300-100-789	60mm Jaw	Long Forward Curved Clamp To be use with a Ø12.5mm Trocar System